





BED ARMCHAIR ACCOMPANYING



- Pediatrics
- Neonatal care
- Maternity
- Hospitalisation
- Aftercare
- Nursing Homes

MARINAGO





MARINAGO - BED ARMCHAIR ACCOMPANYING

STRENGTHS

Safety



Wheel brake leverage

Ergonomics



Retractable castors

Durability




Bed base : suspended electrowelded steel mesh

Comfort



Independent seat from the bed surface
Easy opening mechanism
12 cm mattress
Several bed width available


CHOICE OF COLORS



Granny



Turquoise



Bégonia



Lilas



Souris




Chocolat




Taupe



Nuit



Tomette



Amande



Quetsche




Lavande


35 additional colors available on request



OPTIONS




Cushions for breastfeeding




Thin armrests (-6 cm overall width)


TECHNICAL DETAILS




92




102



60



70



85

h.45

102

Mechanism opening position

Positions

Bed base

Feet

Armrests and backrests

Mattress

Mattress cover

Overall mattress (W x L x D)

Seat

Overall in chair position (W x L x D)

Overall in bed position (W x L x D)

Maximum weight

FR covering M1 cleanable with hospital products, choice of color in color chart

Frontal, bed armchair

2 : seats, extended

Steel Grid, electrowelded connected with springs

4 castors of which 2 retractable

Removable mounted on rails

Hi Density foam 35 Kg /m3 12 cm thickness

FR Polyurethane M1 with steel zipper

Version 60 : 60 x 190 x 12 cm - Version 70 : 70 x 190 x 12 cm

Height 45 cm, width 79 cm, depth 62 cm

Version 60 : 92 x 85 x 102 cm - Version 70 : 102 x 85 x 102 cm

Version 60 : 92 x 85 x 208 cm - Version 70 : 102 x 85 x 208 cm

170 Kg