

Table de gynécologie à hauteur variable / Variable height gynecological table

CARACTÉRISTIQUES STANDARD

Structure acier avec revêtement gris nuage
Hauteur variable électrique de 54 à 93 cm
Dossier réglable par vérin à gaz ou vérin électrique
Coussin de jambes rabattable (suivant modèles)
ou amovible (option)
Paires de poignées
Paire d'étriers ou d'appui-jambes
(suivant modèles)
Support rouleau papier
Tiroir plastique en ABS
Inclinaison manuelle de l'assise de 16°
Patins de réglage sous piètement
Commande par pédale
Capacité patient 200 kg



Capacité / Capacity


200 kg



STANDARD FEATURES

Steel frame with light grey epoxy covering
Electric variable height from 54 to 93 cm
Gas-spring assisted or electric backrest
Foldable leg cushion (according to models)
or removable leg cushion (option)
Pair of hand grips
Pair of stirrups or legrests
(according to models)
Paper roll holder
ABS drawer
Manual seat tilt 16°
Adjustable levellers under base
Pedal control
200 kg patient weight capacity



Références References	Dossier gaz Gas-spring assisted backrest	Dossier électrique Electric backrest	Étriers Stirrups	Appui-jambes Legrests	Sans coussin de jambes Without leg cushion	Coussin de jambes rabattable Foldable leg cushion
BX4110-O1	✓		✓			✓
BX4014-O1	✓			✓	✓	
BX4114-O1	✓			✓		✓
BX4124-O1		✓		✓		✓


Accessoires Accessories

22172



Coussin de jambes amovible
pour BX4014-O1
Removable leg cushion
for BX4014-O1

861



Paire de housses
pour appui-jambes
Pair of covers for legrests





CHOIX DE SELLERIES NON-FEU M1 - PROTECTION ANTIBACTÉRIENNE
CHOICE OF M1 FIRE RETARDANT - ANTIBACTERIAL UPHOLSTERY


Etain
Tin



Marine
Marine



Obscure
Darkness



Raisin
Grape

